

Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test Device (saliva) Package Insert

A RAPID TEST FOR THE QUALITATIVE DETECTION OF NOVEL CORONAVIRUS ANTIGENS IN HUMAN SALIVA.

For professional In Vitro Diagnostic Use Only.

INTENDED USE

The Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test Device (saliva) is an *in vitro* diagnostic test for the qualitative detection Nucleoprotein and spike glycoprotein of Coronavirus Disease 2019 in human Oropharyngeal saliva, using the rapid immunochromatographic method as an aid in the diagnosis of SARS-Cov-2 infections. The identification is based on the monoclonal antibodies specific for the novel coronavirus antigen. It will provide information for clinical doctors to prescribe correct medications.

SUMMARY

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases. Severe acute respiratory syndrome - coronavirus- 2 (SARS-CoV-2) is an enveloped, non-segmented Positive sense rna virus. It is the cause of the Coronavirus-0 disease (COVID-19) common to humans is contagious. SARS-CoV-2 has several structural proteins, including spike (S), envelope (E), membrane (M) and nucleocapsid(N). At present, there are many variants of the Novel coronavirus (SARS-CoV-2), and the N501Y mutation and its approximate variants have attracted attention because their mutation position is located in the spike glycoprotein's receptor-binding domain of the virus, thereby changing the virus infected efficiency. At present, the in-silico analysis shows that the N501Y mutation did not change the primary and tertiary protein structure of the S-RBD protein of the virus, and thus did not change its antigenicity.

PRINCIPLE

The Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test Device (saliva) is an immunochromatographic membrane assay that uses highly sensitive monoclonal antibodies to Novel Coronavirus.

The test strip is composed of the following three parts, namely sample pad, reagent pad and reaction membrane. The reagent membrane contains the colloidal-gold conjugated with the monoclonal antibodies against Novel coronavirus; the reaction membrane contains the secondary antibodies for Novel coronavirus, and the polyclonal antibodies against the mouse globulin, which are pre-immobilized on the membrane.

When the saliva sample is received by the test, the conjugated solution from the reagent pad gets dissolved and migrates along with the saliva. When the Novel coronavirus is present in the saliva sample, a complex is formed between the anti-Novel coronavirus conjugate and the virus will be caught / detected by the specific anti- Novel coronavirus monoclonal coated on the T region. Whether the sample contains the virus or not, the solution continues to migrate to encounter another reagent (an anti-mouse IgG antibody) that binds the remaining conjugates, thereby producing a red line on the region C.

The Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test Device (saliva) product can detect SARS-Cov-2 nucleo-protein (mainly) and the spike glycoprotein of Coronavirus. Thought Elisa methodology, we find our raw material of anti spike protein antibody allows the binding site of 511 to 531 amino acid of the spike protein of SARS-COV-2, which meets the requirements of our original design.

The detection results and sensitivity of the The Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test Device (saliva) are identically accurate when testing for Recombinant spike protein of Novel coronavirus S protein with 319 to 541 amino acid sequence of B.1.1.7(UK variant), B.1.351(SA variant) and standard/wild variants.

The Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test Device (saliva) contains two sets (N+S) of protein detection systems. This allows for a superior antiviral variability detection capability as compared to a single protein spike detection system, which can only detect proteins S or N. This specific double detection design makes The Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test Device (saliva) detection of the nucleo-protein and spike protein of genetic SARS-Cov-2 variants of SARS-CoV-2 much more reliable.

REAGENTS

The reagent membrane contains the colloidal-gold conjugated with the monoclonal antibodies against Novel coronavirus; the reaction membrane contains the secondary antibodies for Novel coronavirus, and the polyclonal antibodies against the mouse globulin, which are pre-immobilized on the membrane.

PRECAUTIONS

- For *in vitro* diagnostic use only.
- Do not use after the expiration date.
- Ensure foil pouch containing test device is not damaged before opening for use.
- Perform test at room temperature 15 to 30°C.
- Wear gloves when hanging the samples, avoid touching the reagent membrane and sample window.
- All samples and used accessories should be treated as infectious and discarded according to local regulations.

• Avoid using bloody samples.

STORAGE AND STABILITY

Store the Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test device (saliva) at room temperature or refrigerated (2-30°C). Do not freeze. All reagents are stable until the expiration dates marked on their outer packaging and buffer vial.

SPECIMEN COLLECTION AND PREPARATION

1. Specimen collection:

The oral fluid specimen should be collected using the collection tools provided with the kit. Follow the detailed Directions for Use below. No other collection tools should be used with this assay. Oral fluid collected at any time of the day may be used.

2. Specimen preparation:

When the saliva is collected, follow the direction to prepare the specimen with buffer provided with the kit.

MATERIALS

Materials provided

- Test device
- Package Insert
- Tube stand*
- Dropper
- Nozzle
- Saliva collect cup/bag
- Extraction buffer
- Extraction tube
- Plastic bag

*The 20-test package contains the tube stand, the 1-test and 5-test package use the test box itself as tube stand.


Materials required but not provided

- Timer


DIRECTIONS FOR USE

Allow the test device, specimen, extraction buffer to equilibrate to room temperature (15-30°C) prior to testing. Do not place anything in the mouth including food, drink, gum, tobacco, water and mouthwash products for at least 10 minutes prior to collection of oral fluid specimen.

1. Spit enough saliva into the saliva collect cup/bag.
2. draw the saliva from the cup with a dropper, transfer 4 drops of saliva to the extraction tube.
3. Take out an extraction tube and a bottle of extraction buffer, remove the extraction buffer bottle cap, add all the extraction buffer into the extraction tube.
4. Take out a nozzle and close into the extraction tube, gently shake the extraction tube vertically for about 5 seconds to allow saliva mix well with extraction buffer.
5. Fold the used cup/bag in half and discard it into the plastic bag as medical waste in accordance with local regulations



6. Remove the test device from the sealed foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed immediately after opening the foil pouch. Put the test device on a clean and flat surface.
7. Transfer 3 drops of sample into the sample well of test device vertically, start the timer.
8. Read the result at 10-20 minutes. Don't interpret the result after 20 minutes.



INTERPRETATION OF RESULTS

(Please refer to the illustration above)

POSITIVE: Two red lines appear. One red line appears in the control region(C), and one red line in the test region(T). The shade of color may vary, but it should be considered positive whenever there is even a faint line.

NEGATIVE: Only one red line appears in the control region(C), and no line in the test region(T). The negative result indicates that there are no Novel coronavirus particles in the sample or the number of viral particles is below the detectable range.

INVALID: No red line appears in the control region(C). The test is invalid even if there is a line on test region(T). Insufficient sample volume or incorrect procedural techniques are the most likely

reasons for control line failure. Review the test procedure and repeat the test using a new test device. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

LIMITATIONS

- The Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test device (saliva) is an acute-phase screening test for qualitative detection. Sample collected may contain antigen concentration below the reagent's sensitivity threshold, so a negative test result does not exclude infection with novel coronavirus.
- The Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test device (saliva) detects viable and non-viable novel coronavirus antigen. Test performance depends on antigen load in the sample and may not correlate with cell culture performed on the same sample. A positive test does not rule out the possibility that other pathogens may be present, therefore, the results must be compared with all other available clinical and laboratory information to make an accurate diagnosis.
- A negative test result may occur if the level of extracted antigen in a specimen is below the sensitivity of the test or if poor quality specimen is obtained.
- Performance of the test has not been established for monitoring antiviral treatment of novel coronavirus.
- Positive test results do not rule out co-infections with other pathogens.
- Negative test results are not intended to rule in other coronavirus infection except the SARS-Cov-2.
- Children tend to shed virus for longer periods of time than adults, which may result in differences in sensitivity between adults and children List.
- The concentration of virus in saliva is greatly affected by factors such as meals, diet, smoking, breath fresheners, etc. Therefore, please strictly follow this manual before collecting samples. A negative result may occur if the concentration of antigen in a specimen is below the detection limit of the test or if the specimen was collected or transported improperly, therefore a negative test result does not eliminate the possibility of SARS-Cov-2 infection, and should be confirmed by viral culture or PCR.

PERFORMANCE CHARACTERISTICS

Clinical Evaluation

Clinical evaluation was performed to compare the results obtained by Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test device (saliva) and PCR test result.

The performance has been determined using 405 Oropharyngeal saliva from individual symptomatic patients with Suspected COVID-19. Oropharyngeal saliva was made refer to the instructions for use of the specimen collection and preparation. All samples were selected and then tested sequentially in a blinded fashion. The performance of the test device was matched with the results of a compared to commercial RT-PCR reagent. The test device showed a sensitivity of 92.9% and a specificity of 99.58%. The results were summarized below.

Table: Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test device (saliva) vs. PCR

| Method | 2019-nCoV Nucleic Acid Test Kit (RT-PCR) | | Total Results | |
|---|--|----------|---------------|----------|
| | Results | Positive | | Negative |
| The Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test device (saliva) | Positive | 157 | 1 | 158 |
| | Negative | 12 | 235 | 247 |
| Total Results | | 169 | 236 | 405 |

Clinical sensitivity = 157/169 = 92.9 % (95%CI*:87.89% to 96.00%)

Clinical specificity = 235/236 = 99.58% (95%CI*:97.39% to >99.99%)

Accuracy: (157+235) / (157+1+12+235) *100% = 96.79% (95%CI* 94.53% to 98.17%)

*Confidence Interval

Limit of Detection (LoD)

| 2019-nCoV Strain Tested | Ready Tech product | | | | |
|---|---|-------------------|----------------------|----------------------|----------|
| Stock 2019-nCoV Concentration | 1 X 10 ⁸ TCID ₅₀ /mL | | | | |
| Dilution | 1/100 | 1/200 | 1/400 | 1/800 | 1/1600 |
| Concentration in Dilution tested (TCID ₅₀ /ml) | 1X10 ⁶ | 5X10 ⁵ | 2.5X 10 ⁵ | 1.25X10 ⁵ | 52.5 |
| Call rates of 20 replicates near cut-off | 100(20/20) | 100(20/20) | 100(20/20) | 95(19/20) | 10(2/20) |
| Limit of detection (LoD) per Virus Strain | 1.25 X 10 ⁵ TCID ₅₀ /mL | | | | |

Cross Reaction

The Cross reaction study conducted By using The Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test Device (saliva) to test sample in which add each concentration of pathogens listed in the following table to the specimens prepared negative samples and 3XLOD positive samples. The results show that the pathogenic listed in the following table have no cross-reaction that would create false positive or negative results for SARS-Cov-2 antigen..

| Virus/Bacteria/Parasite | Strain | Concentration |
|-------------------------|-------------|--|
| MERS-coronavirus | N/A | 72 µg/mL |
| | Type 1 | 1.5 x 10 ⁶ TCID ₅₀ /mL |
| | Type 3 | 7.5 x 10 ⁶ TCID ₅₀ /mL |
| | Type 5 | 4.5 x 10 ⁶ TCID ₅₀ /mL |
| | Type 7 | 1.0 x 10 ⁶ TCID ₅₀ /mL |
| Adenovirus | Type 8 | 1.0 x 10 ⁶ TCID ₅₀ /mL |
| | Type 11 | 2.5 x 10 ⁶ TCID ₅₀ /mL |
| | Type 18 | 2.5 x 10 ⁶ TCID ₅₀ /mL |
| | Type 23 | 6.0 x 10 ⁶ TCID ₅₀ /mL |
| | Type 55 | 1.5 x 10 ⁶ TCID ₅₀ /mL |
| | Influenza A | H1N1 Denver |
| H1N1 WS/33 | | 2.0 x 10 ⁶ TCID ₅₀ /mL |

| | | |
|---|---------------------------------------|--|
| | H1N1 A/Mal/302/54 | 1.5 x 10 ⁸ TCID ₅₀ /mL |
| | H1N1 New Caledonia | 7.6 x 10 ⁸ TCID ₅₀ /mL |
| | H3N2 A/Hong Kong/8/68 | 4.6 x 10 ⁸ TCID ₅₀ /mL |
| Influenza B | Nevada/03/2011 | 1.5 x 10 ⁸ TCID ₅₀ /mL |
| | B/Lee/40 | 8.5 x 10 ⁸ TCID ₅₀ /mL |
| | B/Taiwan/2/62 | 4.0 x 10 ⁸ TCID ₅₀ /mL |
| Respiratory syncytial virus | N/A | 2.5 x 10 ⁸ TCID ₅₀ /mL |
| Legionella pneumophila | Bloomington-2 | 1 x 10 ⁶ PFU/mL |
| | Los Angeles-1 | 1 x 10 ⁶ PFU/mL |
| | 82A3105 | 1 x 10 ⁶ PFU/mL |
| Rhinovirus A16 | N/A | 1.5 x 10 ⁸ TCID ₅₀ /mL |
| Mycobacterium tuberculosis | K | 1 x 10 ⁶ PFU/mL |
| | Erdman | 1 x 10 ⁶ PFU/mL |
| | HN878 | 1 x 10 ⁶ PFU/mL |
| | CDC1551 | 1 x 10 ⁶ PFU/mL |
| | H37Rv | 1 x 10 ⁶ PFU/mL |
| Streptococcus pneumonia | 4752-98 [Maryland (D1)6B-17] | 1 x 10 ⁶ PFU/mL |
| | 178 [Poland 23F-16] | 1 x 10 ⁶ PFU/mL |
| | 262 [CIP 104340] | 1 x 10 ⁶ PFU/mL |
| | Slovakia 14-10 [29055] | 1 x 10 ⁶ PFU/mL |
| Streptococcus pyrogens | Typing strain T1 [NCIB 11841, SF 130] | 1 x 10 ⁶ PFU/ml |
| | Mutant 22 | 1 x 10 ⁶ PFU/ml |
| Mycoplasma pneumoniae | FHstrainofEatonAgent [NCTC 10119] | 1 x 10 ⁶ PFU/ml |
| | 36M129-B7 | 1 x 10 ⁶ PFU/ml |
| | 229E | 1.5 x 10 ⁶ TCID ₅₀ /ml |
| Coronavirus | OC43 | 1.5 x 10 ⁶ TCID ₅₀ /ml |
| | NL63 | 1.5 x 10 ⁶ TCID ₅₀ /ml |
| | HKU1 | 1.5 x 10 ⁶ TCID ₅₀ /ml |
| Human etapneumovirus (hMPV) 3 Type B1 | Peru2-2002 | 1.5 x 10 ⁶ TCID ₅₀ /ml |
| Human Metapneumovirus (hMPV) 16 Type A1 | IA10-2003 | 1.5 x 10 ⁶ TCID ₅₀ /ml |
| Parainfluenza virus | Type 1 | 1.5 x 10 ⁶ TCID ₅₀ /ml |
| | Type 2 | 1.5 x 10 ⁶ TCID ₅₀ /ml |
| | Type 3 | 1.5 x 10 ⁶ TCID ₅₀ /ml |
| | Type 4A | 1.5 x 10 ⁶ TCID ₅₀ /ml |



Interfering Substances Reaction









When tested using the Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test Device (saliva), there was no interference between the device reagents and the Potential interference substances listed in below table that would create false positive or negative results for SARS-Cov-2 antigen.

| Substance | Concentration | Substance | Concentration |
|-----------------------------------|---------------|-------------------------|---------------|
| Mucin | 100µg/mL | Acetylsalicylic acid | 3.0 mM |
| Whole Blood | 5% (v/v) | Ibuprofen | 2.5 mM |
| Biotin | 100µg/mL | Mupirocin | 10 mg/mL |
| Neo-Syneprine (Phenylephrine) | 5%(v/v) | Tobramycin | 10µg/mL |
| Afrin Nasal Spray (Oxymetazoline) | 5%(v/v) | Erythromycin | 50µM |
| Saline Nasal Spray | 5%(v/v) | Ciprofloxacin | 50µM |
| Homeopathic | 5%(v/v) | Ceftriaxone | 110mg/mL |
| Sodium Cromoglycate | 10 mg/mL | Meropenem | 3.7µg/mL |
| Olopatadine Hydrochloride | 10 mg/mL | Tobramycin | 100µg/mL |
| Zanamivir | 5 mg/mL | Histamine Hydrochloride | 100µg/mL |
| Oseltamivir | 10 mg/mL | Peramivir | 1mmol/mL |
| Artemether-lumefantrine | 50µM | Flunisolide | 100µg/mL |
| Doxycycline hyclate | 50µM | Budesonide | 0.64nmol/ L |
| Quinine | 150µM | Fluticasone | 0.3ng/mL |
| Lamivudine | 1 mg/mL | Lopinavir | 3µg/mL |
| Ribavirin | 1 mg/mL | Ritonavir | 3.2mg/mL |
| Daclatasvir | 1 mg/mL | Abidor | 417.8ng/mL |
| Acetaminophen | 150µM | Pooled human mouth wash | N/A |

Dose Hook Effect

Test SARS-CoV-2 wild strain culture medium (concentration 1X10⁷ TCID₅₀/ml) and multiple dilution sample with The Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test Device (saliva), no hook effect was found.

| Symbol | Meaning | Symbol | Meaning |
|--|------------------------------------|---|---------------------------|
|  | In vitro diagnostic medical device |  | Storage temperature limit |

| | | | |
|---|---------------------|---|---|
|  | Manufacturer |  | Authorized representative in the European Community |
|  | Date of Manufacture |  | Use by date |
|  | Do not reuse |  | Consult instruction for use |
|  | Batch code |  | Meet the requirements of EC Directive 98/79/EC |



HANGZHOU REALY TECH CO., LTD.
4th Floor, #12 Building, Eastern Medicine Town,
Xiasha Economic & Technology Development,
310018 Hangzhou, Zhejiang, P. R. China
Website: www.realytech.com



Luxus Lebenswelt GmbH
Kochstr.1.47877, Willich, Germany



Number: 1100003002
Version:2.9
Effective Date:2021-02-09